

ClinicBOXX ClinicBOXX max

Container for
Healthcare waste



Peter Ridley

Waste Systems

ClinicBOXX ClinicBOXX max

Container for Healthcare waste

- UN certified for solid hazardous materials from the following Authorised Packaging Groups: I (X), II (Y), III (Z)
- UN approval no. 1H2W/X.../S/..TGM-...
- Compliant with different hazard levels and UN classifications, e.g. UN 3291
- Approved for transport in accordance with the International Carriage of Dangerous Goods by Road (ADR)

Corporate Identity Options

- Thermal transfer printing
- Hot foiling, e.g. city name, coat of arms, logo, injection moulded texts
- Screen printing

Simple locking system

Simple dual clasp locking system for easy manual opening and closing



Certified bins for bio-medical waste

e.g. Clinical Waste in accordance with UN 3291 and other Hazardous Medical Waste categories. Wastes classified as non-hazardous e.g. AHP's – local regulations apply

Patented Lid

Patented lid design with two locking clasps and lid seal, allows quick and easy opening and closing by the operative. The container's large aperture allows filling without spillages.




Easy to handle

Ergonomic, asymmetric handle with a 30mm diameter for optimum handling with gloves.



Chip nest

Integrated chip nest for transponder chips under the comb bar, the container offers excellent possibilities of identification and registration of bin data.



Heavy Duty Wheels

ClinicBOXX® and ClinicBOXX®max are extra mobile thanks to its vulcanised rubber wheels and heavy duty hubs, with option for the wider BigFootWheel.



Quick-Release Wheels

Empty ClinicBOXX® and ClinicBOXX®max containers can be nested inside each other (without wheels), as wheels are fitted with a quick release mechanism. This saves a lot of space for storage and transport.



HENKEL


In addition, all current regulations for the collection of Hazardous Waste Materials must be observed by the operator and user.

Quality:

ClinicBOXX and ClinicBOXXmax corresponds with standard EN 840. Manufactured under Quality Management System ISO 9001:2015 and compliant with the stringent automotive standard IATF 16949:2016. Regular quality monitoring is undertaken by our in-house laboratory and independent EU Test Institutes.



	Volume in litres (Nominal)	Height(h) Width(w) Depth(d) mm	Wheels Ø mm	Weight kg (approx.)
<p>ClinicBOXX</p>	120 litres	960/ 498/ 554	200	11,5
<p>ClinicBOXXmax</p>	200 litres	1045/ 580/ 674	200	15



Peter Ridley Waste Systems
 Unit , Eastlands Industrial Estate,
 Eastlands Road, Leiston,
 Suffolk, IP16 4LL
 Tel: 01728 660755
 Email: office@peterridley.co.uk
 Web: www.peterridley.co.uk



FUTURE SOLUTIONS TODAY.